

AOHUA



VETERINARY

Endoscopic Accessories

Munich, Germany

Tel: +49(0)81 04/89908-10

Address: Fichtenstrasse 27, Brunnthal, Munich, Germany



Aohua Endoscopy Co.,Ltd.

Location: No.66, Lane133,Guangzhong Road, Minhang District, Shanghai, P.R.China

Postcode: 201108

Email: marketing@aohua.com

www.aohua.com/vet.cn



ISO 9001/ISO 13485 Certified



CONTENTS

MAKE ENDOSCOPY

DIAGNOSIS AND TREATMENT

MORE ACCESSIBLE AND CONVENIENT

CONTENTS

01

Biopsy Forceps

01

02

Grasping Forceps

03

03

Accessories for
HF Surgery

06

04

Cleaning Brush

07

05

Veterinary Stent

08

06

Other Accessories

09

07

Accessories Kit

10



01 Biopsy Forceps



Biopsy Forceps, Sterile, Pack of 10

Disposable

	Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
	Oval Cup without Spike Uncoated	JRQ-Y-1816	φ1.8x1600	≥2.0
		JRQ-Y-1823	φ1.8x2300	≥2.0
		JRQ-Y-2323	φ2.3x2300	≥2.8
		JRQ-Y-2340	φ2.3x4000	≥2.8
	Oval Cup with Spike Uncoated	JRQ-Z-2323	φ2.3x2300	≥2.8
		JRQ-Z-2340	φ2.3x4000	≥2.8
	Aligator Cup Uncoated	JRQ-Y-1816-EA	φ1.8x1600	≥2.0
		JRQ-Y-1823-EA	φ1.8x2300	≥2.0
		JRQ-Y-2323-EA	φ2.3x2300	≥2.8
		JRQ-Y-2340-EA	φ2.3x4000	≥2.8
	Aligator Cup with Spike	JRQ-Z-2323-EA	φ2.3x2300	≥2.8
		JRQ-Z-2340-EA	φ2.3x4000	≥2.8

Unsterile, Pack of 1

Reusable

	Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
	Long Oval Cup not Fenestrated	AH100-10-160	φ1.0x1600	≥1.2
	Long Oval Cup Serrated	AH-118-18-160	φ1.8x1600	≥2.0
		AH-118-18-230	φ1.8x2300	≥2.0
		AH-118-23-230	φ2.3x2300	≥2.8
		AH-118-23-400	φ2.3x4000	≥2.8
		AH-118-34-400	φ3.4x4000	≥3.7
	Long Oval Cup Fenestrated	AH-100-18-160	φ1.8x1600	≥2.0
		AH-100-18-230	φ1.8x2300	≥2.0
		AH-100-23-230	φ2.3x2300	≥2.8
		AH-100-23-400	φ2.3x4000	≥2.8
		AH-100-30-400	φ3.0x4000	≥3.7
	Long Oval Cup Fenestrated with Spike	AH-101-18-160	φ1.8x1600	≥2.0
		AH-101-18-230	φ1.8x2300	≥2.0
		AH-101-23-230	φ2.3x2300	≥2.8
		AH-101-23-400	φ2.3x4000	≥2.8
		AH-101-30-400	φ3.0x4000	≥3.7
	Oval Cup with fine teeth	AH-136-18-160	φ1.8x1600	≥2.0
		AH-136-18-230	φ1.8x2300	≥2.0
		AH-136-23-230	φ2.3x2300	≥2.8
		AH-136-23-400	φ2.3x4000	≥2.8
	Round Cup fenestrated	AH-106-18-160	φ1.8x1600	≥2.0
		AH-106-18-230	φ1.8x2300	≥2.0
		AH-106-23-230	φ2.3x2300	≥2.8
		AH-106-23-400	φ2.3x4000	≥2.8
		AH-106-30-400	φ3.0x4000	≥3.7

Cytology brush, Sterile, Pack of 10

Disposable

	Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
	Straight With sheath	JRS-A φ1.1×1200-PN	φ1.1x1200	≥1.2
		JRS-A φ1.8×1200-PN	φ1.8x1200	≥2.0
		JRS-A φ1.8×2300-PN	φ1.8x2300	≥2.0
		JRS-A φ1.8×4000-PN	φ1.8x4000	≥2.0

Grasping Forceps



Grasping Forceps, Sterile, Pack of 10

Disposable

	Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	with opening width(mm)	Working Channel (mm)
	Rat teeth and Alligator type	JRY-QF-1821	φ1.8x2100	≥8.75	≥2.0
	Alligator type	JRY-QE-2321	φ2.3x2100	≥7	≥2.8
		JRY-QE-2340	φ2.3x4000	≥7	≥2.8
	Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)		Working Channel (mm)
	Rat teeth type	JRY-QD-2321	φ2.3x2100	≥7	≥2.8
		AH-109-10-160	φ1.0x1600		≥1.2
		AH-109-18-160	φ1.8x1600		≥2.0
	Alligator Forceps	JRY-QD-2340	φ2.3x4000	≥7	≥2.8
		AH-109-18-230	φ1.8x2300		≥2.0
		AH-109-23-230	φ2.3x2300		≥2.8
		JRY-QH-2321	φ2.3x2100	≥14.4	≥2.8
		AH-109-30-400	φ3.0x4000		≥3.7
	Pelican Type	JRY-QH-2340	φ2.3x4000	≥14.4	≥2.8
		AH-108-10-160	φ1.0x1600		≥2.0
		AH-108-18-160	φ1.8x1600		≥2.0
	Rat Tooth, 2x1 teeth	AH-108-18-230	φ1.8x2300		≥2.8
		AH-108-23-230	φ2.3x2300		≥2.8
		AH-108-23-400	φ2.3x4000		≥3.7
		AH-108-30-400	φ3.0x4000		≥2.0

Stone Basket, Sterile, Pack of 5

Disposable

Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
 Diamond shape	JRY-WD-1818-25	φ1.8x1800	≥2.0
	JRY-WD-1823-25	φ1.8x2300	≥2.0
	JRY-WD-2423-35	φ2.4x2300	≥2.8
	JRY-WD-2440-35	φ2.4x4000	≥2.8

Grasping Forceps, Unsterile, Pack of 1

Reusable

Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
 Alligator Forceps	AH-109-10-160	φ1.0x1600	≥1.2
	AH-109-18-160	φ1.8x1600	≥2.0
	AH-109-18-230	φ1.8x2300	≥2.0
	AH-109-23-230	φ2.3x2300	≥2.8
	AH-109-23-400	φ2.3x4000	≥2.8
	AH-109-30-400	φ3.0x4000	≥3.7
 Rat Tooth 2x1 teeth	AH-108-10-160	φ1.0x1600	≥1.2
	AH-108-18-160	φ1.8x1600	≥2.0
	AH-108-18-230	φ1.8x2300	≥2.0
	AH-108-23-230	φ2.3x2300	≥2.8
	AH-108-23-400	φ2.3x4000	≥2.8
	AH-108-30-400	φ3.0x4000	≥3.7
 3 prongs Grasper	AH-112-10-160	φ1.0x1600	≥1.2
	AH-112-18-160	φ1.8x1600	≥2.0
	AH-112-18-230	φ1.8x2300	≥2.0
	AH-112-23-230	φ2.3x2300	≥2.8
	AH-112-23-400	φ2.3x4000	≥2.8
 3 prongs Retriever round	AH-114-18-160	φ1.8x1600	≥2.0
	AH-114-18-230	φ1.8x2300	≥2.0
	AH-114-23-230	φ2.3x2300	≥2.8
	AH-114-23-400	φ2.3x4000	≥2.8
 4 prongs Retriever sharp	AH-115-18-160	φ1.8x1600	≥2.0
	AH-115-18-230	φ1.8x2300	≥2.0
	AH-115-23-230	φ2.3x2300	≥2.8
	AH-115-23-400	φ2.3x4000	≥2.8
 Foregin Body Retriever 2x2 teeth	AH-110-23-230	φ2.3x2300	≥2.8
	AH-110-23-400	φ2.3x4000	≥2.8



Forceps, Unsterile, Pack of 1

Reusable

Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
 Foregin Body Retriever 2x1 teeth	AH-111-18-160	φ1.8x1600	≥2.0
	AH-111-18-230	φ1.8x2300	≥2.0
	AH-111-23-230	φ2.3x2300	≥2.8
	AH-111-23-400	φ2.3x4000	≥2.8
 Grasping Forceps Serrated	AH-113-18-160	φ1.8x1600	≥2.0
	AH-113-18-230	φ1.8x2300	≥2.0
	AH-113-23-230	φ2.3x2300	≥2.8
	AH-113-23-400	φ2.3x4000	≥2.8
 Grasping Forceps for nutrition probes	AH-131-18-160	φ1.8x1600	≥2.0
	AH-131-18-230	φ1.8x2300	≥2.0
	AH-131-25-230	φ2.5x2300	≥2.8
	AH-131-25-400	φ2.5x4000	≥2.8
 Polycatch – Forceps	AH-132-25-230	φ2.5x2300	≥2.8
 Grasping Forceps Rat Tooth 2:1 teeth cross serrated	AH-148-23-230	φ2.3x2300 Opening of jaws: 14,9 mm Length of jaws: 8,4 mm	≥2.8
	AH-148-23-400	φ2.3x4000 Opening of jaws: 14,9 mm Length of jaws: 8,4 mm	≥2.8

Accessories for HFE Surgery



Snare, Sterile, Pack of 5

Disposable

Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	loop width(mm)	Working Channel (mm)
	JRD-QA-1823-25	$\phi 1.8 \times 2300$	25	≥ 2.0
	JRD-QA-2423-35	$\phi 2.4 \times 2300$	35	≥ 2.8
	JRD-QA-2440-35	$\phi 2.4 \times 4000$	35	≥ 2.8

Injection Needle, Sterile, Pack of 5

Disposable


Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)	Needle's diameter	Needle's length (mm)
	JRZ-1816-2504	$\phi 1.8 \times 1600$	≥ 2.0	25G	4
	JRZ-1823-2504	$\phi 1.8 \times 2300$	≥ 2.0	25G	4
	JRZ-2323-2305	$\phi 2.3 \times 2300$	≥ 2.8	23G	5
	JRZ-2340-2305	$\phi 2.3 \times 4000$	≥ 2.8	23G	5

Aspiration catheters, Sterile, Pack of 10

Disposable

Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Working Channel (mm)
	ACT120S	$\phi 1.0 \times 1200$	≥ 1.2
	ACT200S	$\phi 1.8 \times 2000$	≥ 2.0
	ACT230S	$\phi 2.3 \times 2300$	≥ 2.8
	ACT400S	$\phi 2.3 \times 4000$	≥ 2.8

Cleaning Brush



Cleaning Brushes, Sterile, Pack of 5

Disposable


Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Bristle Diameter (mm)	Working Channel (mm)
 Double end Cleaning Brushes	JRX1012-0210	φ1.0x1200	1.5x10	≥1.2
	JRX1812-0320	φ1.8x1200	3x12	≥2.0
	JRX1823-0320	φ1.8x2300	3x12	≥2.0
	JRX1823-0520	φ1.8x2300	5x12	≥2.8
	JRX1823-0520	φ1.8x4000	5x12	≥2.8
	JRX1840-0620	φ1.8x4000	6x20	≥3.7

Cleaning Brushes, Sterile, Pack of 1

Reusable

Type	Model NO.	Max inserton Porton Diameter &Working Length(mm)	Bristle Diameter (mm) x Length (mm)	Working Channel (mm)
 Double end Cleaning Brushes	JRF1823-0512	φ1.8x2300	5x12	≥2.8
	JRF1840-0512	φ1.8x4000	5x12	≥2.8
 Valves Brushes	JRF95-1040	φ1.0x95	10x40	/

Veterinary Stent



Veterinary Tracheal Stent (straight shape), Niti-alloy Wire With Introducer

Disposable



Model NO.	Diameter(mm)	Working Channel (mm)
GS-8-11-08030	8	30
GS-8-11-08050	8	50
GS-8-11-08070	8	70
GS-8-11-08080	8	80
GS-8-11-08090	8	90
GS-8-11-10030	10	30
GS-8-11-10050	10	50
GS-8-11-10060	10	60
GS-8-11-10070	10	70
GS-8-11-10080	10	80
GS-8-11-10090	10	90
GS-8-11-12040	12	40
GS-8-11-12050	12	50
GS-8-11-12060	12	60
GS-8-11-12070	12	70
GS-8-11-12090	12	90
GS-8-11-14030	14	30
GS-8-11-14040	14	40
GS-8-11-14060	14	60
GS-8-11-14070	14	70
GS-8-11-14080	14	80
GS-8-11-14100	14	100

Other Accessories

Mucus extractors, individually packed and Sterile, Pack of 5

Disposable



Type

BAL / Tracheal wash catheters can become blocked with thick mucus, collection traps ensure you retrieve the best sample to aid in your diagnosis. Single chamber mucus extractors are;

Gastroscopy Hood

Reusable



Type

In order to reduce the risk of damage to the Cardiac sphincter, Esophagus and Pharynx during foreign body retrieval our protective hood can be fitted to the majority of scopes to ensure safe extraction of a foreign body via gastroscopy. Hood Measurements: 40mm at Distal end, 8.3mm at proximal end, length of 75mm.

Guttural Pouch Probe

Reusable



Diameter(mm)

Working Channel (mm)

1.8

1800

2.4

1800

2.4

2250

Accessories Kit



Disposable

Reusable

VET-3512



QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
☐	Cytology Brush	Straight With sheath	JRS-A $\Phi 1.1 \times 1200$ -PN	1.1	1200
☐	Aspiration Catheters		ACT120S	1.0	1200
☐	Biopsy Forceps	Long Oval Cup	AH100-10-160	1.0	1600
☐	Grasping Forceps	Alligator	AH-109-10-160	1.0	1600
☐	Grasping Forceps	Rat Tooth 2x1 teeth	AH-108-10-160	1.0	1600
☐	Grasping Forceps	3 prongs Grasper	AH-112-10-160	1.0	1600

VET-6010



QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
☐	Biopsy Forceps	Oval Cup	JRQ-Y-1816	1.8	1600
☐	Biopsy Forceps	Aligator Cup	JRQ-Y-1816-EA	1.8	1600
☐	Cytology Brush		JRS-A $\Phi 1.8 \times 1200$ -PN	1.8	1200
☐	Aspiration Catheters		ACT200S	1.8	2000
☐	Grasping Forceps	Rat teeth & Aligator type	JRY-QF-1821	1.8	2100
☐	Stone basket	Diamond shape	JRY-WD-1818-25	1.8	1800
☐	Snare	Oval	JRD-QA-1823-25	1.8	2300
☐	Injection Needle		JRZ-1816-2504	1.8	1600

Accessories Kit

Disposable

Reusable

VET-8015

	QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
	☐	Biopsy Forceps	Oval Cup	JRQ-Y-1823	1.8	2300
	☐	Biopsy Forceps	Aligator Cup	JRQ-Y-1823-EA	1.8	2300
	☐	Cytology Brush		JRS-A Φ1.8×2300-PN	1.8	2300
	☐	Grasping Forceps	Rat teeth & Aligator type	JRY-QF-1821	1.8	2100
	☐	Stone basket	Diamond shape	JRY-WD-1823-25	1.8	2300
	☐	Snare	Oval	JRD-QA-1823-25	1.8	2300
	☐	Injection Needle		JRZ-1823-2504	1.8	2300

VET-8028

	QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
	☐	Biopsy Forceps	Oval Cup	JRQ-Y-2323	2.3	2300
	☐	Biopsy Forceps	Oval Cup with Spike	JRQ-Z-2323	2.3	2300
	☐	Biopsy Forceps	Aligator Cup	JRQ-Y-2323-EA	2.3	2300
	☐	Biopsy Forceps	Aligator Cup with Spike	JRQ-Z-2323-EA	2.3	2300
	☐	Cytology Brush		JRS-A Φ1.8×2300-PN	1.8	2300
	☐	Grasping Forceps	Rat teeth	JRY-QD-2321	2.3	2100
	☐	Grasping Forceps	Aligator type	JRY-QE-2321	2.3	2100
	☐	Grasping Forceps	Pelican type	JRY-QT-2321	2.3	2100
	☐	Stone basket	Diamond shape	JRY-WD-2423-35	2.4	2300
	☐	Snare	Oval	JRD-QA-2423-35	2.4	2300
	☐	Injection Needle		JRZ-2323-2305	2.3	2300

Disposable

Reusable

VET-9215

	QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
	☐	Biopsy Forceps	Oval Cup	JRQ-Y-2323	2.3	2300
	☐	Biopsy Forceps	Oval Cup with Spike	JRQ-Z-2323	2.3	2300
	☐	Biopsy Forceps	Aligator Cup	JRQ-Y-2323-EA	2.3	2300
	☐	Biopsy Forceps	Aligator Cup with Spike	JRQ-Z-2323-EA	2.3	2300
	☐	Cytology Brush		JRS-A Φ1.8×2300-PN	1.8	2300
	☐	Grasping Forceps	Rat teeth	JRY-QD-2321	2.3	2100
	☐	Grasping Forceps	Aligator type	JRY-QE-2321	2.3	2100
	☐	Grasping Forceps	Pelican type	JRY-QT-2321	2.3	2100
	☐	Stone basket	Diamond shape	JRY-WD-2423-35	2.4	2300
	☐	Snare	Oval	JRD-QA-2423-35	2.4	2300
	☐	Injection Needle		JRZ-2323-2305	2.3	2300
	☐	Aspiration Catheters		ACT230S	2.3	2300

VET-9833/VET-1330/VET-1335

	QTY	Type	Features	Model No.	Diameter (mm)	Working Length(mm)
	☐	Biopsy Forceps	Oval Cup	JRQ-Y-2340	2.3	4000
	☐	Biopsy Forceps	Oval Cup with Spike	JRQ-Z-2340	2.3	4000
	☐	Biopsy Forceps	Aligator Cup	JRQ-Y-2340-EA	2.3	4000
	☐	Biopsy Forceps	Aligator Cup with Spike	JRQ-Z-2340-EA	2.3	4000
	☐	Cytology Brush		JRS-A Φ1.8×4000-PN	1.8	4000
	☐	Grasping Forceps	Rat teeth	JRY-QD-2340	2.3	4000
	☐	Grasping Forceps	Aligator type	JRY-QE-2340	2.3	4000
	☐	Grasping Forceps	Pelican type	JRY-QT-2340	2.3	4000
	☐	Stone basket	Diamond shape	JRY-WD-2440-35	2.4	4000
	☐	Snare	Oval	JRD-QA-2440-35	2.4	4000
	☐	Injection Needle		JRZ-2340-2305	2.3	4000
	☐	Aspiration Catheters		ACT400S	2.3	4000